

RS-BLUE

0.05% TRYPAN BLUE

0.05% TRYPAN BLUE, 1.5 ML

TISSUE STAINING

	RS-BLUE
Ingredients	<ul style="list-style-type: none"> ▶ trypan blue ▶ sodium chloride ▶ sodium phosphate monobasic ▶ sodium phosphate dibasic ▶ water for injectables
Dye Concentration	0.045 - 0.055 %w/V
Hydrogenionic activity (pH)	7.00 - 7.40 pH units
Osmolality	255 - 345 mOsm/kg
Osmolarity	254 - 344 mOsm/L
Bacterial Endotoxins	< 0.5 EU/ml
Cytotoxicity Test	Non-cytotoxic according to ISO 10993-5, direct contact test performed on two cell lines: - retinal pigment epithelial cells (ARPE-19) - mouse fibroblasts cells (BALB 3T3 clone A31)
Form	sterile liquid
Volume	1.5 ml
Packaging	vial
Sterilization	moist heat according to EN 556-1 and EN ISO 17665-1
Storage temperature	+15°C / +32°C
Expiry date	36 months
Intended use	reversible tissue staining during ophthalmic surgery
Class	medical device
GMDN code	P 45180
CND code	Q 02030201
MD registration number	3608/R
Medical device class	Ila
CE mark	CE 0477

RS-BLUE is 0.05% trypan blue, 1.5 ml, packaged in vial, for reversible staining of human tissues during surgery. Product warrants a prompt and effective staining of the capsule during cataract surgery, thus allowing a clear visualization of tissues without any risks for patients.

Ref.: RSB 001-00

Contents:

- ▶ 6 RS-Blue vials, 1.5 ml
- ▶ 1 package insert
- ▶ 12 labels for product traceability

GS1-128 barcode

ALCHIMIA

Viale Austria, 14
35020 Ponte S.Nicolò (PD)
Italy

VAT IT 00063370282

Phone +39 049 8962074 / 64
Fax +39 049 8962071

www.alchimiasrl.com
info@alchimiasrl.com